

Shun Guider



Introducer Sets



Femoral Sheath

Code	Size(Fr)	Length(cm)	Needle Gauge(Ga)	Guidewire Diameter(inch)	Guidewire Length(cm)
Standard Kit					
615401	4F	9	20	0.021	45
* 615501	5F	11	18	0.035	45
* 615601	6F	11	18	0.038	45
* 615701	7F	11	18	0.038	45
* 615801	8F	11	18	0.038	45
615901	9F	11	18	0.038	45
615101	10F	11	18	0.038	45
615908	6F	11	18	0.035	45
With hydrophilic coating					
615972	4F	9	20	0.021	45
615973	5F	11	18	0.035	45
615974	6F	11	18	0.038	45
615975	7F	11	18	0.038	45
615976	8F	11	18	0.038	45
615977	9F	11	18	0.038	45
615978	10F	11	18	0.038	45

Long introducer sheath

Code	Size(Fr)	Length(cm)	Needle Gauge(Ga)	Guidewire Diameter(inch)	Guidewire Length(cm)
Standard Kit					
615502	5F	23	18	0.035	70
615602	6F	23	18	0.038	70
615702	7F	23	18	0.038	70
615802	8F	23	18	0.038	70
615902	9F	23	18	0.038	70
615102	10F	23	18	0.038	70
615911	6F	45	18	-	-
615912	7F	45	18	-	-
615913	8F	45	18	-	-
With hydrophilic coating					
615984	5F	23	18	0.035	70
615980	6F	23	18	0.038	70
615981	7F	23	18	0.038	70
615982	8F	23	18	0.038	70
615983	9F	23	18	0.038	70
615985	10F	23	18	0.038	70
615914	6F	45	18	-	-
615915	7F	45	18	-	-
615916	8F	45	18	-	-

Radial Sheath

Code	Size(Fr)	Length(cm)	Needle Gauge(Ga)	Guidewire Diameter(inch)	Guidewire Length(cm)
Standard Kit					
615419	4F	7	20*3.8	0.021	45
615420	4F	11	20*3.8	0.021	45
* 615423	5F	7	20*3.8	0.021	45
* 615540	5F	11	20*3.8	0.021	45
615968	4F	7	21*3.8	0.018	45
615969	4F	11	21*3.8	0.018	45
615963	6F	7	21*3.8	0.018	45
615964	5F	7	21*3.8	0.018	45
* 615542	6F	7	20*3.8	0.021	45
* 615543	6F	11	20*3.8	0.021	45
615966	6F	11	21*3.8	0.018	45
615967	5F	11	21*3.8	0.018	45
615996	5F	11	20*3.2	0.021	45
615821	6F	11	20*3.2	0.021	45
615545	7F	7	20*3.8	0.021	45
615546	7F	11	20*3.8	0.021	45
615970	7F	7	21*3.8	0.018	45
615971	7F	11	21*3.8	0.018	45
With hydrophilic coating					
615988	4F	7	21*3.8	0.018	45
615989	4F	11	21*3.8	0.018	45
615990	6F	7	21*3.8	0.018	45
615991	6F	11	21*3.8	0.018	45
615992	5F	7	21*3.8	0.018	45
615993	5F	11	21*3.8	0.018	45
615994	7F	7	21*3.8	0.018	45
615995	7F	11	21*3.8	0.018	45
615840	4F	7	20*3.8	0.021	45
615841	4F	11	20*3.8	0.021	45
* 615842	6F	7	20*3.8	0.021	45
* 615843	6F	11	20*3.8	0.021	45
* 615844	5F	7	20*3.8	0.021	45
* 615845	5F	11	20*3.8	0.021	45
615846	7F	7	20*3.8	0.021	45
615847	7F	11	20*3.8	0.021	45

Note: The code with * mark means main selling items