



Intravascular Lithotripsy Generator and Connector Cable



Product features

- Difficult to open the heavy calcification
- High risk to dissection
- Safe and efficiency

Ordering information

Code	Model	Length of IVL generator(mm)	Width of IVL generator(mm)	Height of IVL generator(mm)
629968	IVLG-I	292	77	200

Peripheral Intravascular Lithotripsy Catheter



Ordering information

Code	Model (PIVL14)	Balloon nominal diameter (mm)	Balloon length (mm)	Balloon treatment pressure (atm)	Nominal pressure (atm)	Rated bursting pressure (atm)	Compatible introducer sheath (Fr)
629964	PIVL14-25400	2.50	40	4	6	10	6
631143	PIVL14-25600	2.50	60	4	6	10	6
629965	PIVL14-30400	3.00	40	4	6	10	6
631144	PIVL14-30600	3.00	60	4	6	10	6
629966	PIVL14-35400	3.50	40	4	6	10	6
629956	PIVL14-35600	3.50	60	4	6	10	6

Ordering information

Code	Model (PIVL14)	Balloon nominal diameter (mm)	Balloon length (mm)	Balloon treatment pressure (atm)	Nominal pressure (atm)	Rated bursting pressure (atm)	Compatible introducer sheath (Fr)
631145	PIVL14-35800	3.50	80	4	6	10	6
629967	PIVL14-40400	4.00	40	4	6	10	6
629957	PIVL14-40600	4.00	60	4	6	10	6
631146	PIVL14-40800	4.00	80	4	6	10	6
629958	PIVL14-45600	4.50	60	4	6	10	6
631147	PIVL14-45800	4.50	80	4	6	10	6
629959	PIVL14-50600	5.00	60	4	6	10	6
631148	PIVL14-50800	5.00	80	4	6	10	6
629960	PIVL14-55600	5.50	60	4	6	10	7
631149	PIVL14-55800	5.50	80	4	6	10	7
629961	PIVL14-60600	6.00	60	4	6	10	7
631150	PIVL14-60800	6.00	80	4	6	10	7
629962	PIVL14-65600	6.50	60	4	6	10	7
631151	PIVL14-65800	6.50	80	4	6	10	7
629963	PIVL14-70600	7.00	60	4	6	10	7
631152	PIVL14-70800	7.00	80	4	6	10	7
631153	PIVL14-80600	8.00	60	4	6	10	7
631154	PIVL14-80800	8.00	80	4	6	10	7

*Compatible guide wire: 0.014 inch

Ordering information


Code	Model (PIVL18)	Balloon nominal diameter (mm)	Balloon length (mm)	Balloon treatment pressure (atm)	Nominal pressure (atm)	Rated bursting pressure (atm)	Compatible introducer sheath (Fr)
631155	PIVL18-25400	2.50	40	4	6	10	6
631156	PIVL18-25600	2.50	60	4	6	10	6
631157	PIVL18-30400	3.00	40	4	6	10	6
631158	PIVL18-30600	3.00	60	4	6	10	6
631159	PIVL18-35400	3.50	40	4	6	10	6
631160	PIVL18-35600	3.50	60	4	6	10	6
631170	PIVL18-35800	3.50	80	4	6	10	6
631171	PIVL18-40400	4.00	40	4	6	10	7
631172	PIVL18-40600	4.00	60	4	6	10	7
631173	PIVL18-40800	4.00	80	4	6	10	7
631174	PIVL18-45600	4.50	60	4	6	10	7

Ordering information

Code	Model (PIVL18)	Balloon nominal diameter (mm)	Balloon length (mm)	Balloon treatment pressure (atm)	Nominal pressure (atm)	Rated bursting pressure (atm)	Compatible introducer sheath (Fr)
631175	PIVL18-45800	4.50	80	4	6	10	7
631176	PIVL18-50600	5.00	60	4	6	10	7
631177	PIVL18-50800	5.00	80	4	6	10	7
631178	PIVL18-55600	5.50	60	4	6	10	7
631179	PIVL18-55800	5.50	80	4	6	10	7
631180	PIVL18-60600	6.00	60	4	6	10	7
631181	PIVL18-60800	6.00	80	4	6	10	7
631182	PIVL18-65600	6.50	60	4	6	10	7
631183	PIVL18-65800	6.50	80	4	6	10	7
631184	PIVL18-70600	7.00	60	4	6	10	7
631185	PIVL18-70800	7.00	80	4	6	10	7
631186	PIVL18-80600	8.00	60	4	6	10	7
631187	PIVL18-80800	8.00	80	4	6	10	7

*Compatible guide wire: 0.018 inch

Coronary Intravascular Lithotripsy Catheter



Ordering information

Code	Model	Balloon nominal diameter (mm)	Balloon length (mm)	Balloon treatment pressure (atm)	Nominal pressure (atm)	Rated bursting pressure (atm)	Compatible introducer sheath (Fr)
629952	CIVL25012	2.50	12	4	6	10	6
630584	CIVL27512	2.75	12	4	6	10	6
629953	CIVL30012	3.00	12	4	6	10	6
630585	CIVL32512	3.25	12	4	6	10	6
629954	CIVL35012	3.50	12	4	6	10	6
629955	CIVL40012	4.00	12	4	6	10	6

Ordering information

Code	Model	Balloon nominal diameter (mm)	Balloon length (mm)	Balloon treatment pressure (atm)	Nominal pressure (atm)	Rated bursting pressure (atm)	Compatible introducer sheath (Fr)
630560	CIVL25015	2.50	15	4	6	10	6
630586	CIVL27515	2.75	15	4	6	10	6
630561	CIVL30015	3.00	15	4	6	10	6
630587	CIVL32515	3.25	15	4	6	10	6
630562	CIVL35015	3.50	15	4	6	10	6
630563	CIVL40015	4.00	15	4	6	10	6



SHUNMEI MEDICAL Co., Ltd.

Head Office : R301 and R401, Building B, No. 8, Jinlong First Road, Baolong Industrial Zone, Longgang District, Shenzhen, 518116 Guangdong, P.R. China

Factory: Yifa 3rd road, Yifa Industrial Zone, Pingtan Town, Huiyang District, Huizhou, Guangdong, China

TEL: 0086-0752-3306949

Email: nora@shunmed.com

Web: www.shunmed.com

