


CTO Coronary Balloon Dilatation Catheter



* Note: Balloon catheters with diameters of 1.0mm and 1.5mm only have a radiopaque mark in the middle of the balloon.

Product feature

- CTO balloon superior pass-through design for CTO lesions;
- Special balloon capsule process design, resulting in a smaller passing diameter;
- Enhanced push bar design, bringing better push and support;
- The soft and smooth tip enhances the passability and minimizes the damage to the vascular wall;
- Can be repeated multiple expansion, slender side to enhance passability;
- Balloon diameter increases with the increase of pressure, stable performance and bursting pressure value.

Clinical application department

- Department of Cardiology
- Chest pain center
- Intervention department
- Stroke center

Diameter(mm)	Length(mm)						
	6	8	9	10	12	15	16
1.0	630343	630535	630344	630536	630386	630537	/
1.25	/	630538	630539	630540	630541	630542	/
1.5	/	630543	630544	630545	630546	630547	630548
1.75	/	630549	630550	630551	630552	630553	630554