


Shunmei[®]

Medical technology creates future


Shun Ski Trails

Hydrophilic Coated Guidewire


Straight tip



Angled tip



J tip



1.5mm/3.0mm
Diameter



PRODUCT CODE	WIRE TYPE	DIAMETER	FLEXIBLE TIP LENGTH	TIP SHAPE	LENGTH
626355	Standard	0.032in	3cm	Straight tip	150cm
626356					180cm
626357					260cm
626440				Angled tip	150cm
626441					180cm
626442					260cm
626455				J tip (3mm)	150cm
626456					180cm
626457					260cm
632070				J tip(1.5mm)	150cm
632071					180cm
632072					260cm
626358		0.035in		Straight tip	150cm
626359					180cm
626360					260cm
626443				Angled tip	150cm
626444					180cm
626445					260cm
626458				J tip (3mm)	150cm
626459					180cm
626460					260cm
632073				J tip(1.5mm)	150cm
632074					180cm
632075					260cm
626361		0.038in		Straight tip	150cm
626362					180cm
626363					260cm
626446				Angled tip	150cm
626447					180cm
626448					260cm
626461				J tip (3mm)	150cm
626462					180cm
626464	260cm				
632076	J tip(1.5mm)		150cm		
632077			180cm		
632078			260cm		

PRODUCT CODE	WIRE TYPE	DIAMETER	FLEXIBLE TIP LENGTH	TIP SHAPE	LENGTH	
632034	Stiff	0.032in	3cm	Straight tip	150cm	
632035					180cm	
632036					260cm	
632037				Angled tip	150cm	
632038					180cm	
632039					260cm	
632040				J tip (3mm)	150cm	
632041					180cm	
632042					260cm	
632043				J tip(1.5mm)	150cm	
632044					180cm	
632045					260cm	
632046				Straight tip	150cm	
632047					180cm	
632048					260cm	
632049		Angled tip		150cm		
632050				180cm		
632051				260cm		
632052		0.035in		3cm	J tip (3mm)	150cm
632053						180cm
632054						260cm
632055		J tip(1.5mm)		150cm		
632056				180cm		
632057				260cm		
632058		Straight tip		150cm		
632059				180cm		
632060				260cm		
632061		Angled tip		150cm		
632062				180cm		
632063				260cm		
632064		0.038in		3cm	J tip(3mm)	150cm
632065						180cm
632066						260cm
632067	J tip(1.5mm)	150cm				
632068		180cm				
632069		260cm				

 **SHUNMEI MEDICAL CO.,LTD**

Add: R401 of Building B, No.8 1st jinlong Road, Baolong Industrial Zone, Longgang District, Shenzhen City, Guangdong China.

Branch Office: Yifa Industrial Zone, Dushi Village, Pingtan Town, Huiyang District, Huizhou, China.

Mobile/Whatsapp: 86-183-4436 9007

86-135-1076 5957

Email: info@shunmed.com

Website: www.shunmed.com

