

Slender™

Guide
Extension
Catheter



Shunmei®

SHUNMEI AND ITS SUBSIDIARIES

2026-03-25

SHUNMEI MEDICAL CO.,LTD

Add: No.8 jinlong 1st street, Baolong Industrial Zone,
Longgang District, Shenzhen, China.

Branch Office: Yifa Industrial Zone, Dushi Village,
Pingtan Town, Huiyang District, Huizhou, China.

Mobile/Whatsapp: 86-183-4436 9007 +86 183 4435 9973

E-mail: nora@shunmed.com nora1@shunmed.com

Website: www.shunmed.com



APP

SHUNMEI SP. Z O.O.

SHUNMEİ MEDİKAL VE TIBBİ CİHAZLAR
TİCARET ANONİM ŞİRKETİ



Shunmei

Search



Shunmei®

Slender™



**External Markings +
Dual Radiopaque Rings**
Easy Advancement Assessment



Optional Side-Hole Design,
Less Coronary Ischemia,
Better Tolerance



**Braided +
Polymer Structure**
Excellent Support



OD: **5F–8F**
Length: **25mm/40mm**



Tapered Distal Shaft,
Seamless Catheter
Segment Transition



PTFE Hydrophobic
Coating On Proximal Shaft,
Smooth Delivery

Guide Extension Catheter Models and Sizes

Code	Model	Effetive Length(cm)	Cather Length(cm)	Adapted guiding catheter inner diameter(mm)	Cather Inner Diameter(mm)	Catheter outer diameter(mm)	Number of side holes/aper ture
635746	6F-25	150	25	ID≥1.78	1.45	1.71±0.06	0
635747	6F-40	150	40	ID≥1.78	1.45	1.71±0.06	0
635748	7F-25	150	25	ID≥1.98	1.6	1.86±0. 1	0
635749	8F-25	150	25	ID≥2.24	1.83	2. 11±0. 1	0
635750	6F-25 SH8	150	25	ID≥1.78	1.45	1.71±0.06	8/0.6mm
635751	6F-40 SH8	150	40	ID≥1.78	1.45	1.71±0.06	8/0.6mm
635752	7F-25 SH8	150	25	ID≥1.98	1.6	1.86±0. 1	8/0.6mm
635753	8F-25 SH8	150	25	ID≥2.24	1.83	2. 11±0. 1	8/0.6mm

Intended Use :

- Support:**

Provide additional guidance support to the lesion site;

- Extension:**

Prolonging the operational range of therapeutic devices;

- Delivery:**

Offering a smooth pathway for specialized devices (such as drug-coated balloons, stents, etc);